

Biotechnology Explorer™

CAPTIVATING SCIENCE EDUCATION

BIO-RAD

Professional Development

ELISA Immuno Explorer™: Antibodies in Agriculture From Mad Cow to GMOs



ELISA Immuno Explor[™]

Instructors



Stan Hitomi

Coordinator – Math & Science
San Ramon Valley Unified School District
Danville, CA

Kirk Brown

Lead Instructor, Edward Teller Education Center
Science Chair, Tracy High School
and Delta College, Tracy, CA

Sherri Andrews, Ph.D.

Curriculum and Training Specialist
Bio-Rad Laboratories

Essy Levy, M.Sc.

Curriculum and Training Specialist
Bio-Rad Laboratories

Why Teach ELISA?



- **Hands-on Immunology**
- **Tangible results**
- **Laboratory extensions**
- **Real-world connections**
- **Link to careers and industry**
- **Standards-based:**
 - **One lesson integrates multiple standards**
 - **Health sciences**
 - **Immunology**
 - **Biodefense**
 - **Immune response** – antibody/antigen interactions
 - **Disease** – infection, detection, transmission

Scientific Inquiry

- Tapping nature's tool kit to solve human problems
- Use of immunodetection to hunt for proteins
- Use of positive and negative experimental controls
- Interpretation of experimental results

Genetics

- DNA > RNA > protein — antibodies
- Antibody structure and function
- Antibody production via genetic recombination

Cell and Molecular Biology

- Immune response
- Manufacturing antibodies
- Virology and immunology

Chemistry of Life

- Enzyme-substrate interactions
- Protein structure and function
- Properties of antigens and antibodies

Evolution

- Animal immune systems response
- HIV mutation and evolution
- Viral drug resistance
- Biowarfare in nature

Environmental and Health Science

- HIV, mad cow disease, and bird flu testing
- Epidemiology and biodefense
- Drug, pregnancy, and GMOs testing
- Soil, water, air testing

ELISA Immuno Explorer Kit Advantages

- **Lab completed in a 45 min period**
- **Supplies for 48 students (12 workstations)**
- **Comprehensive and flexible curriculum**
- **Compelling real-world links**
- **Striking results**
- **Cost effective**
- **Classroom Safe**



Workshop Time Line

- **Introduction**
- **Antigen Detection by ELISA**
- **Ways the ELISA-Immuno Explorer Kit can be used**
- **Real-World Examples**

ELISA

Enzyme-Linked Immuno**s**orbant Assay




- **Mammalian immune system**
- **Antibody specificity**
- **Biology’s “magic bullet”**
- **Evolved over millions of years**
- **Harness nature’s tool kit**
- **Imagine the applications!**

Links to the Real World




- **Mad Cow Disease, SARS, HIV**
- **GMO**
- **Drug and steroid testing**
- **Pregnancy / Reproduction**
- **Biodefense**
- **Cancer treatment**

Immune Response





ELISA Antibody Structure



ELISA ANIMATION

ELISA Procedures Overview



ELISA Kit Workstation Inventory



Reagents:

Yellow tubes	Test samples	2
Violet tube (+)	Positive control	1
Blue tube (-)	Negative control	1
Green tube (PA)	Primary antibody	1
Orange tube (SA)	Secondary antibody	1

Lab Equipment and Supplies:

Microplate strips, pipettor, pipette tips, transfer pipette, wash buffer, paper towels, marking pen

Laboratory Quick Guide

Laboratory Quick Guide

Antigen Detection ELISA

Student Workstation Checklist

One workstation serves 4 students.

Item (Label)	Contents	Number	(✓)
Yellow tubes	Student test samples (0.25 ml)	4	<input type="checkbox"/>
Violet tube (+)	Positive control (0.5 ml)	1	<input type="checkbox"/>
Blue tube (-)	Negative control (0.5 ml)	1	<input type="checkbox"/>
Green tube (PA)	Primary antibody (1.5 ml)	1	<input type="checkbox"/>
Orange tube (SA)	Secondary antibody (1.5 ml)	1	<input type="checkbox"/>
Brown tube (SUB)	Enzyme substrate (1.5 ml)	1	<input type="checkbox"/>
12-well microplate strips		2	<input type="checkbox"/>
50 µl fixed-volume micropipet or 20–200 µl adjustable micropipet		1	<input type="checkbox"/>
Yellow tips		10–20	<input type="checkbox"/>
Disposable plastic transfer pipet		1	<input type="checkbox"/>
70–80 ml wash buffer in beaker	Phosphate buffered saline with 0.05% Tween 20	1	<input type="checkbox"/>
Large stack of paper towels		2	<input type="checkbox"/>
Black marking pen		1	<input type="checkbox"/>

1. Label the yellow tubes with each student's initials.
2. Label your 12-well strip. On each strip label the first 3 wells with a "+" for the positive controls and the next 3 wells with a "-" for the negative controls. Label the remaining wells with your and your partner's initials (3 wells each).



3. Use a **fresh** pipet tip to transfer 50 µl of the positive control (+) into the three "+" wells.
4. Use a **fresh** pipet tip to transfer 50 µl of the negative control (-) into the three "-" wells.
5. Transfer 50 µl of each of your team's samples into the appropriately initialed three wells, using a **fresh** pipet tip for each sample.
6. Wait 5 minutes while the proteins in the samples bind to the plastic wells.



PROTOCOL II
Antigen Detection ELISA


Step One

Label and add controls


- Obtain a test-sample
- Label the 12-well strip:
 - First 3 wells: positive controls “+”
 - Next 3 wells: negative controls “-”
 - Remaining wells to identify test-samples



- Add 50 ul of positive control to 1st 3 wells
- Add 50 ul of negative control to 2nd 3 wells
- Add 50ul of the student samples to the appropriately labeled wells
- Wait 5 minutes for the antigen to bind



Microplate Strips




- **Microplate strips are made of polystyrene**
- **Hydrophobic side chains in amino acids bind to the polystyrene wells**



- **No coating is needed**

Step Two

WASH




- **Remove samples from wells by firmly tapping them on a paper towel**
- **Discard the top paper towel**
- **Using a disposable transfer pipette wash wells with wash buffer**
- **Remove wash buffer by firmly tapping the wells on a paper towel**
- **Discard the top paper towel**
- **Repeat wash step**


Step Three

Add (PA) Primary Antibody

- Add 50 ul of the primary antibody (PA) to all 12 wells
- Samples are left in wells for 5 minutes
- After 5 minutes WASH 2X




Wash Buffer




- **Wash buffer contains phosphate buffer saline (PBS) to keep antibodies in a stable environment that helps keep their structure**
- **Also contains Tween 20: a nonionic detergent removes non-specifically bound proteins and coats wells that acts as a blocking agent to reduce background**
- **Antibody will only bind to the antigen**

Step Four

Wash antibody and add enzyme-linked secondary antibody (SA)







- Wash the primary antibody from polystyrene wells as before
- WASH 2X
- Add 50ul of the enzyme-linked secondary antibody to each well
- Wait 5 minutes




Antibody Specificity

- **Secondary antibody (enzyme-linked antibody) will only bind to the primary antibody (serum antibody)**
- **Secondary antibody specifically recognizes the constant region of the primary antibody**
- **In which wells do you predict this is happening?**




Step Five


Add enzyme substrate (SUB)



- **Wash the enzyme-linked secondary antibody from polystyrene wells as before**
- **Using a disposable transfer pipette wash wells with wash buffer**
- **WASH 3X**
- **Add 50ul of the enzyme substrate to each well**
- **Wait 5 minutes**
- **positive samples will begin to turn blue**



What are the reagents?





Antigen: Chicken gamma globulin

Primary antibody (PA): Polyclonal anti-chicken antibody made by rabbits

Secondary antibody (enzyme-linked) SA: Polyclonal anti-rabbit antibody made by goats linked (conjugated) to horseradish peroxidase (HRP)

Enzyme substrate (SUB): 3,3',5,5' – tetramethylbenzidine (TMB) – a colorless solution that when oxidized by HRP turns blue

ELISA Kit Results




Ways The ELISA Kit Can Be Used

Protocol	Type of ELISA	Real-World Application
I	Tracking outbreaks of disease	HIV, SARS, smallpox & anthrax
II	Detecting antigens	GMO, BSE, pregnancy, drugs, (and all the above)
III	Detecting antibodies in serum	HIV, Lyme disease, smallpox and West Nile virus

ELISA test for Transmissible Spongiform Encephalopathies (TSEs)


- Uses differences in diseased prions vs. normal prions to prepare sample.
- Proteinase K only digests normal, **not diseased**, prions .
- ELISA tests for any prion protein

Prion Proteins (PrP^{res} and PrP^{sens})



PrP^{sens}

- Proteinase K **sensitive**
- Soluble in detergent



PrP^{res}

- Proteinase K **resistant**
- Aggregates in detergent

TSE test sample preparation

1. **Sample brain tissue**
2. **Homogenize brain tissue**
3. **Digest with Proteinase K**
(normal prions are digested, diseased prions are resistant)
4. **Concentrate**
5. **Denature Proteinase K**
6. **Perform ELISA**



Protocol II: Antigen Detection ELISA

Protocol - ELISA on simulated animal brain samples



Real-World Application – TSE Test

Tube Description	Actual Tube Contents	Simulated Tube Contents
Student samples	Antigen or PBS	Processed brain
Primary antibody	Primary antibody	Antibody against prion protein
Secondary antibody	Secondary antibody	HRP-linked antibody against primary antibody
Positive control	Antigen	Synthesized peptide with prion sequence
Negative control	PBS	Buffer

Real-world Applications of Antibodies

Applications

- Dipstick tests/ELISA
- Immunostaining
- Western blotting

Agricultural Uses

- Crop-specific disease diagnosis
- Animal disease diagnosis
- Detection of GM crops
- Basic research



ELISA to test for GMOs



“Genetically Modified Organism (GMO)”

an organism in which the genetic material has been altered in a way that does not occur naturally by mating and/or natural recombination

DNA → **RNA** → **Protein**

- **ELISA can help farmers separate their GMO grain lots from non-GMO grain lots.**
- **ELISA tests are used to identify specific proteins**
 - **Delta-endotoxin Cry1Ab from Bt11**
 - **glyphosate from Round-up (RR)**

How to test for GMOs



ELISA:

Test for presence of proteins expressed from genetic modifications

Pro: Quick, inexpensive, low tech

Con: Crop specific, protein stability


PCR:

Test for presence of inserted foreign DNA


Pro: ID different GM crops, DNA stability

Con: Expensive, timely

Example: Pregnancy Test



Positive



Negative

