

(S)2-Propanol WARNING/ATTENTION: LIQUID AND VAPOR Causes eye irritation. Hazard if swallowed - may irritate if swallowed or vomited. Prolonged or repeated skin contact may cause drying, cracking, or irritation. High vapor concentrations may cause dizziness and inhibition of the eyes or respiratory tract. Keep away from heat, sparks and flame. Take precautionary measures against static discharge. Avoid breathing high vapor concentrations. Avoid contact with eyes and prolonged skin contact. Use only with adequate ventilation. Keep container closed. Wash thoroughly after handling. SEE MATERIAL SAFETY DATA SHEET.

(C)2-Propanol ATTENTION/AVVERTIMENTO: LIQUIDE ET VAPEUR IRRITANTS. Provoque une irritation des yeux. Nocif en cas d'ingestion ou de contact prolongé. Le contact répété ou prolongé avec la peau peut provoquer un séchage, un fendillement et une irritation. Des concentrations élevées de vapeurs peuvent provoquer de la confusion et de la nausée. Prenez des mesures de précaution contre les décharges électrostatiques. Gardez le récipient fermé. Évitez le contact avec la peau. Utilisez uniquement avec une ventilation adéquate. Gardez le récipient fermé. Lavez-vous soigneusement après manipulation. VOIR FICHE SIGNALÉTIQUE.



Avantor Performance Materials, Inc., USA

(S)2-Propanol DANGER. Highly flammable liquid and vapor. Causes serious eye irritation. May cause respiratory irritation. May cause dizziness or disorientation. Keep away from heat/sparks/open flame/hot surfaces. No smoking. Wear protective gloves/protective clothing/face protection/foot protection. Call a POISON CENTER or doctor/physician if you feel unwell. Keep container tightly closed.

(R)2-Propanol. Perigo. Líquido e vapor altamente inflamáveis. Pode ser nocivo em caso de ingestão e por prolongado. Causa irritação ocular séria. Pode causar irritação respiratória. Pode causar tonturas e orientação. Mantenha afastado de calor/faixas/chamas. Não fume. Use níveis de proteção adequados de proteção/botão protetor de proteção facial. Chame o CENTRO DE ASSOCIAÇÃO TOXICOLÓGICA ou um médico se não se sentir bem. SE NECESSÁRIO, procure atendimento médico imediatamente. Não provoque vômito. Manter a embalagem bem fechada.

(C)2-Propanol. 危險。高度可燃性液體和蒸氣。吞嚥時進入呼吸道可傷害眼睛。造成嚴重眼刺激。可能引起呼吸系統刺激。可能引起頭暈或定向障礙。遠離火源、火花、明火、熱表面。嚴禁吸菸。穿戴適當的手套/防護服/防護面罩/防護鞋。如感不適，請向毒物中心諮詢。如感不適，立即向醫務中心或醫生。不得誘吐。保持容裝物密封。

(T)2-Propanol. 危險。高度可燃性液體和蒸氣。吞嚥時進入呼吸道可傷害眼睛。造成嚴重眼刺激。可能引起呼吸系統刺激。可能引起頭暈或定向障礙。遠離火源、火花、明火、熱表面。嚴禁吸菸。穿戴適當的手套/防護服/防護面罩/防護鞋。如感不適，請向毒物中心諮詢。如感不適，立即向醫務中心或醫生。不得誘吐。保持容裝物密封。



4L
Isopropyl Alcohol
Suitable for Liquid Chromatography and UV Spectrophotometry

ChromAR®
CH₃CHOHCH₃
FW 60.10
Batch No.
Manufactured Date: YYYY/MM/DD
Expiration Date: YYYY/MM/DD
Meets ACS Specifications
For Laboratory Research or Manufacturing Use
Product Information (not certification) Appearance (clear, colorless liquid)
For a Certificate of Lot Analysis, visit www.avanormaterials.com

UN1219

FLASH POINT: 12°C (54°F)
(Closed Cup)
DOT Name: ISOPROPANOL

CAS NO: 67-63-0

Made in USA

WARNING: This product contains chemicals known to the state of California to cause cancer.

MPI #
Rev #

E24

AWANTOR
PERFORMANCE MATERIALS
Avantor Performance Materials, Inc.
3477 Corporate Parkway, Suite 200
Center Valley, PA 15204
(610) 573-2600
www.avanormaterials.com



Lot/Batch Number

Shelf Life / Expiration / Use by Date

Supplier Product Number